



PRODUCT SPECIFICATION SHEET

Haemophilus Test Agar Base (DM1273)

Intended Use

Haemophilus Test Agar Base (DM1273) is recommended for the susceptibility testing of *Haemophilus influenzae*.

Product Summary and Explanation

Haemophilus species are nutritionally fastidious in nature requiring either exogenous hemin (X-factor) or NAD (Vfactor) or both.⁽¹⁾ *Haemophilus influenzae* require complex media for growth. These complex media have aggravated the routine susceptibility testing of *H. influenzae* because of antagonism between some essential nutrients and certain antimicrobial agents. In 1966, Bauer, Kirby and others developed a standardized procedure for the antimicrobial susceptibility testing of common, rapidly growing bacteria in which Mueller Hinton Agar was selected as the test medium. This medium is not satisfactory for fastidious organisms such as some streptococci, gonococci and *Haemophilus* species.⁽²⁻⁴⁾ Also, addition of blood to Mueller Hinton Agar to supply the essential growth nutrients makes the medium opaque, rendering it unsuitable for antimicrobial susceptibility testing.

Haemophilus Test Agar Base, studied by Jorgensen et al^(5, 6) has been specifically formulated for the susceptibility testing of *Haemophilus influenzae*. This medium has similar composition as Mueller Hinton Agar, with the addition of yeast extract and added growth supplements. Haemophilus Test Agar Base is simple, transparent and poses minimum risk of antagonism of antimicrobial agents.⁽⁵⁾ Haemophilus Test Agar Base is also recommended by the United States National Committee for Clinical Laboratory Standards (NCCLS) for both dilution and disc diffusion assays.⁽⁷⁾ A major advantage of HTM Agar compared with Mueller Hinton Chocolate Agar is optical clarity, permitting zone diameter measurements from the bottom of the dish as is the standard test procedure for non-fastidious organisms on Mueller Hinton Agar. Comparisons with Mueller-Hinton Chocolate Agar have shown an overall agreement of 99.6%.⁽⁸⁾ Furthermore, HTM Agar contains low levels of thymidine and is, therefore, suitable for testing trimethoprim/sulfamethoxazole.

Principles of the Procedure

Haemophilus Test Agar Base contains beef infusion and casein acid hydrolysate, which provides carbon, nitrogen and other essential nutrients required for growth of organisms. Yeast extract serves as a source of B complex vitamins. Starch acts as a protective colloid against toxic substances present in the medium.

Formula / Liter

Ingredients	Gms / Liter
Beef infusion from	300.00
Casein acid hydrolysate	17.50
Yeast extract	5.00
Starch	1.50
Agar	17.00
Final pH: 7.4 ± 0.2 at 25°C	
Formula may be adjusted and/or supplemented as required to meet performance specifications	

Precautions

1. For Laboratory Use only.
2. IRRITANT. Irritating to eyes, respiratory system, and skin.

Directions

1. Suspend 21.5 grams of the medium in 500 ml of distilled water.
2. Heat to boiling, to dissolve the medium completely.
3. Autoclave at 121°C, 15 psi pressure, for 15 minutes / validated cycle.
4. Cool to 50°C and aseptically add the rehydrated contents of 1 vial of Haemophilus Growth Supplement (MS085).
5. Mix well and pour into sterile Petri plates.





PRODUCT SPECIFICATION SHEET

Quality Control Specifications

Dehydrated Appearance	Cream to yellow homogeneous free flowing powder
Prepared Medium	Light amber coloured, clear to slightly opalescent gel forms in Petri plates
Reaction of 4.3% solution	pH 7.4 ± 0.2 at 25°C
Gel Strength	Firm, comparable with 1.7% Agar gel

Expected Cultural Response: Cultural characteristics observed with added Haemophilus Growth Supplement (MS085) in 5-7% carbon dioxide after an incubation at 35-37°C for 18-24 hours.

Sr. No.	Organisms	Results to be achieved		
		Inoculum (CFU)	Growth	Recovery
1.	<i>Haemophilus influenzae</i> ATCC 49766	50-100	good-luxuriant	≥70%
2.	<i>Enterococcus faecalis</i> ATCC 29212	50-100	good-luxuriant	≥70%
3.	<i>Streptococcus pyogenes</i> ATCC 19615	50-100	good-luxuriant	≥70%
4.	<i>Neisseria meningitidis</i> ATCC 13090	50-100	good-luxuriant	≥70%
5.	<i>Staphylococcus aureus</i> ATCC 25923	50-100	good-luxuriant	≥70%

The organisms listed are the minimum that should be used for quality control testing.

Test Procedure

1. Prepare a Gram stain before starting susceptibility testing to confirm culture purity and to confirm tentative identification of *Haemophilus*.
2. Use several well-isolated colonies taken directly from an overnight (preferably 20-24 hours) Chocolate Agar plate (DM314) as the source of the inoculum.
3. A rapid β-lactamase test should be utilized for rapid detection of strains that are resistant to penicillin, ampicillin or amoxicillin.
4. Prepare a suspension of the test organism in Mueller Hinton Broth (DM173) or 0.9% saline. This suspension should be adjusted to the turbidity of the 0.5 McFarland standard using a photometric device. This suspension will contain approximately $1-4 \times 10^8$ CFU/mL. Care must be exercised in preparing this suspension because higher inoculum concentrations may lead to false-resistant results with some β-lactam antibiotics, particularly when β-lactamase-producing strains of *H. influenzae* are tested.
5. Consult the product literature or the CLSI Approved Standard M2⁽⁹⁾ for details on plate inoculation and use of antimicrobial discs.

Bauer-Kirby procedure:

1. This procedure is based on the diffusion through an agar gel of antimicrobial substances which are impregnated on paper discs.
2. Earlier methods made use of discs from high and low antimicrobial concentrations and the presence or absence of inhibition zones for their interpretation, in this method discs with single concentration of antimicrobial agent and zone diameters are correlated with minimal inhibitory concentrations (MICs).
3. In the test procedure, a standardized suspension of the organism is swabbed over the entire surface of the medium. Paper discs impregnated with specified amounts of antibiotic or other antimicrobial agents are then placed on the surface of the medium, the plate is incubated and zones of the inhibition around each disc are measured.
4. The determination as to whether the organism is susceptible, resistant or intermediate in its response to the agent is made by comparing zone diameters obtained to those provided with CLSI document M2.⁽⁹⁾
5. Various factors have been identified as influencing disc diffusion susceptibility tests. These include the medium, agar depth, disc potency, inoculum concentration, pH and beta-lactamase production by test organisms.⁽⁹⁻¹²⁾

Results

1. Examine the plates after 16-18 hours of incubation. A confluent "lawn" of growth should be obtained. If only isolated colonies grow, the inoculum was too light and the test should be repeated.
2. Measure the diameter of the zones of complete inhibition (as judged by the unaided eye), including the diameter of the disc, to the nearest whole millimeter, using calipers, a ruler, or a template prepared for this purpose. The measuring device is held on the back of the plate, which is held over a black, non-reflecting background and illuminated from above. The endpoint should be taken as the area showing no obvious visible growth that can be





PRODUCT SPECIFICATION SHEET

detected with the unaided eye. Disregard faint growth of tiny colonies which can be detected with difficulty near the edge of the obvious zone of the inhibition.

3. Consult the product literature or the CLSI Approved Standard M2⁽⁹⁾ for details on interpretation of results.

Storage

Store the sealed bottle containing the dehydrated medium at 2 - 30°C. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light.

Expiration

Refer to the expiration date stamped on the container. The dehydrated medium should be discarded if not free flowing, or if the appearance has changed from the original color. Expiry applies to medium in its intact container when stored as directed.

Limitations of the Procedure

1. For identification, organisms must be in pure culture. Morphological, biochemical and/or serological tests should be performed for final identification.
2. Consult appropriate texts for detailed information and recommended procedures.

Packaging

Product Name : Haemophilus Test Agar Base

Product Code : DM1273

Available Pack sizes : 500gm

References

1. Murray P. R., Baron J. H., Tenover F. C. and Tenover F. C., (Eds.), 2003, *Manual of Clinical Microbiology*, 8th Ed., American Society for Microbiology, Washington, D.C.
2. Bauer A. W., Kirby W. M., Sherris J. C. and Tenover F. C., 1966, *Am. J. Clin. Pathol.* 45:493.
3. Ryan K. J., Schoenknecht F. D., and Kirby W. M., 1970, *Hospital Practice*, 5:91.
4. Barry A. L., Garcia F., and Thrupp L. D., 1970, *Am. J. Clin. Pathol.*, 53 :149.
5. Jorgensen J. H., Redding J. S., Maher L. A. and Howell A. W., 1987, *J. Clin. Microbiol.*, 25:2105.
6. Jorgensen J. H., Howell A. W., and Maher L. A., *J. Clin. Microbiol.*, 28:985.
7. NCCLS Documents M2-A4 Vol. 10. No. 7. and M7-A2 Vol. 10. No 8.
8. Evans G., Marsik F., Thompson L. and Fowler J. (1990) Abstracts of ASM Meeting 1990 C-252.
9. Clinical and Laboratory Standards Institute. 2006. Approved standard: M2-A9. Performance standards for antimicrobial disk susceptibility tests, 9th ed. CLSI, Wayne, Pa.
10. Neumann, Sahn, Thornsberry and McGowan. 1991. *Cumitech 6A*, New developments in antimicrobial agent susceptibility testing, A practical guide. Coord. ed., McGowan. American Society for Microbiology, Washington, D.C.
11. Ericsson and Sherris. 1971. *Acta. Pathol. Microbiol. Scand. Sect. B, Suppl.* 217:1.
12. Jorgensen, Turnidge and Washington. 1999. *In Murray, Baron, Tenover and Tenover (ed.), Manual of clinical microbiology*, 7th ed. American Society for Microbiology, Washington, D.C.

Further Information

For further information please contact your local MICROMASTER Representative.



MICROMASTER LABORATORIES PRIVATE LIMITED

Unit 38/39, Kalpataru Industrial Estate,

Off G.B. Road, Near 'R-Mall', Thane (W) - 400607. M.S. INDIA.

Ph: +91-9320126789/98333630009/9819991103

Email: sales@micromasterlab.com

DM1273PI, Rev.0, 01.08.2008



PRODUCT SPECIFICATION SHEET



Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

