



# PRODUCT SPECIFICATION SHEET

## Clostridium Difficile Agar Base (DM062)

### Intended Use

Clostridium Difficile Agar Base (DM062) is recommended for selective isolation of *Clostridium difficile* from faeces.

### Product Summary and Explanation

*Clostridium difficile* is a pathogenic *Clostridium* affecting the bowel. The spectrum of disease caused by *Clostridium difficile* ranges from pseudomembranous colitis (PMC) through antibiotic associated colitis (AAC). It is now one of the most commonly detected enteric pathogens, and an important cause of nosocomial infections in hospitals and nursing homes.<sup>(1)</sup> It also includes chronic inflammatory bowel diseases, post-operative diarrhoea and non-antibiotic associated diarrhoea.<sup>(2)</sup> The organism has been isolated from diverse natural habitats, including soil, hay, sand, dung from various large mammals, and from dog, cat, rodent and human faeces. Smith and King<sup>(3)</sup> first reported the presence of *C. difficile* in human infections. In 1979, George et al. developed a medium called CCFA (cycloserine-cefoxitin-fructose agar), based on the Egg Yolk Agar formula of McClung and Toabe with fructose replacing glucose.<sup>(4)</sup> Clostridium Difficile Agar is a modification of the original CCFA formulation.<sup>(1)</sup> The medium was made inhibitory to the accompanying flora by the addition of the selective agents namely, D-cycloserine and cefoxitin. This medium does not contain neutral red indicator, as in the original formulation, as it is recommended for use with sheep or horse blood.<sup>(4)</sup> Clostridium Difficile Agar Base is used for the primary isolation of *C. difficile* from faecal specimens.

### Principles of the Procedure

Clostridium Difficile Agar Base contains proteose peptone which provides nitrogen, vitamins, and amino acids in Clostridium Difficile Agar. Fructose is the fermentable carbohydrate used to enhance recovery and growth of *C. difficile*. The Phosphates are buffering agents in this medium. Magnesium Sulfate is a source of inorganic ions to stimulate growth. Sodium Chloride maintains the osmotic balance of the medium. The selective agents D-cycloserine and cefoxitin used in the medium inhibit the growth of majority of *Enterobacteriaceae* and also *Enterococcus faecalis*, Staphylococci, gram-negative anaerobic bacilli and *Clostridium* species other than *C. difficile*, which may be found abundantly in faecal samples. Addition of 7% v/v horse blood to the base increases the recovery of *C. difficile* and also increases its colony size.

### Formula / Liter

Ingredients	Gms / Liter
Proteose peptone	40.00
Disodium phosphate	5.00
Monopotassium phosphate	1.00
Magnesium sulphate	0.10
Sodium chloride	2.00
Fructose	6.00
Agar	15.00
Final pH: 7.4 ± 0.2 at 25°C	
Formula may be adjusted and/or supplemented as required to meet performance specifications	

### Precautions

1. For Laboratory Use only.





## PRODUCT SPECIFICATION SHEET

2. IRRITANT. Irritating to eyes, respiratory system, and skin.

### Directions

1. Suspend 34.55 grams of the medium in 500ml of distilled water.
2. Heat gently to boiling to dissolve the medium completely.
3. Autoclave at 121°C, 15 psi pressure, for 15 minutes / validated cycle.
4. Cool to 50°C. Aseptically add rehydrated contents of 1 vial of Clostridium Difficile Supplement (MS011) together with 7% (v/v) defibrinated Horse blood or Sheep blood.
5. Mix well and pour into sterile Petri plates.

### Quality Control Specifications

Dehydrated Appearance	Cream to yellow homogeneous free flowing powder
Prepared Medium	Basal medium: Light amber coloured clear to slightly opalescent gel After addition 7% v/v defibrinated horse blood: Cherry red coloured, opaque gel forms in Petri plates
Reaction of 6.91% solution	pH 7.4 ± 0.2 at 25°C
Gel Strength	Firm, comparable with 1.5% Agar gel

**Expected Cultural Response:** Cultural characteristics observed under anaerobic condition with added Clostridium Difficile Supplement (MS011) and 7% v/v defibrinated horse blood, after an incubation at 35-37°C for 48 hours.

Sr. No.	Organisms	Results to be achieved			
		Inoculum (CFU)	Growth	Recovery	Colour of Colony
1.	<i>Clostridium difficile</i> ATCC 11204	50-100	good-luxuriant	>=50%	greyish-white
2.	<i>Shigella flexneri</i> ATCC 12022	>=10 <sup>3</sup>	inhibited	0%	
3.	<i>Escherichia coli</i> ATCC 25922	>=10 <sup>3</sup>	inhibited	0%	
4.	<i>Staphylococcus aureus</i> ATCC 25923	>=10 <sup>3</sup>	inhibited	0%	

The organisms listed are the minimum that should be used for quality control testing.

### Test Procedure

1. Spread a part of the faecal sample on the medium to obtain isolated colonies. Incubate the plates anaerobically at 37°C for 18 - 48 hours.
2. Refer to appropriate references for a complete discussion on the isolation and identification of *C. difficile* and other anaerobic bacteria.<sup>(1,5)</sup>

### Results

1. *C. difficile* forms grayish white, irregular, raised and opaque colonies, 4-6 mm in diameter, after 48 hours.
2. Typical gram stain morphology of *C. difficile* may not be seen in colonies taken from this medium due to the presence of antibiotics.
3. Subculture on Blood Agar to obtain characteristic morphology.
4. *C. difficile* colonies will not exhibit the typical fluorescence and colour of colony on this medium whereas other Clostridia can give fluorescence. Therefore, for complete identification and confirmation, other tests like gram staining, morphology, biochemicals, specific cytotoxin and clinical observation should be carried out.





## PRODUCT SPECIFICATION SHEET

---

5. Refer to appropriate references and standard test procedures for interpretation of results.

### Results

Examine plates for fungal colonies exhibiting typical color and morphology. Biochemical tests and serological procedures should be performed to confirm findings.

### Storage

Store the sealed bottle containing the dehydrated medium at 10 - 30°C. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light.

### Expiration

Refer to the expiration date stamped on the container. The dehydrated medium should be discarded if not free flowing, or if the appearance has changed from the original color. Expiry applies to medium in its intact container when stored as directed.

### Limitations of the Procedure

1. Clostridium Difficile Agar does not contain Neutral Red indicator because it is designed for use with horse blood.
2. For identification, organisms must be in pure culture. Morphological, biochemical and/or serological tests should be performed for final identification.
3. Consult appropriate texts for detailed information and recommended procedures.

### Packaging

**Product Name : Clostridium Difficile Agar Base**

**Product Code : DM062**

**Available Pack sizes : 500gm**

### References

1. Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (eds.). 1995. Manual of clinical microbiology, 6<sup>th</sup> ed. American Society for Microbiology, Washington, D.C.
2. Collee J. G., Fraser A. G., Marmion B. P., Simmons A., (Eds.), Mackie and McCartney, Practical Medical Microbiology, 14th Ed., Churchill Livingstone.
3. Smith L. D. S. and King E. O., 1962, J. Bacteriol., 84:65.
4. George W. L., Sutter V. L., Citron D., and Finegold S. M., 1979, J.Clin. Microbiol., 9:214.
5. Isenberg, H. D. (ed.). 1992. Clinical microbiology procedures handbook. American Society for Microbiology, Washington, D.C.

### Further Information

For further information please contact your local MICROMASTER Representative.



**MICROMASTER LABORATORIES PRIVATE LIMITED**

Unit 38/39, Kalpataru Industrial Estate,

Off G.B. Road, Near 'R-Mall', Thane (W) - 400607. M.S. INDIA.

Ph: +91-22-25895505, 4760, 4681. Cell: 9320126789.

Email: [micromaster@micromasterlab.com](mailto:micromaster@micromasterlab.com)




[sales@micromasterlab.com](mailto:sales@micromasterlab.com)

DM062PSS,QAD/FR/024,Rev.00/01.01.2018





## PRODUCT SPECIFICATION SHEET

	Checked By	Approved By
 01.01.2018	 01.01.2018	 01.01.2018
Microbiologist	Head Quality Control	Head Quality Assurance

### Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

