



PRODUCT SPECIFICATION SHEET

Acetamide Agar (Twin Pack) (DM444)

Intended Use

Acetamide Agar (DM444) is recommended for differentiation of *Pseudomonas aeruginosa* and other non-fermentative gram-negative bacteria.

Product Summary and Explanation

Acetamide Agar is formulated as per the recommendation of Standard Methods for the Examination of Water and Wastewater.⁽¹⁾ Assimilation studies by Gilardi and others using basal mineral media showed that acetamide was utilized by a wide variety of nonfermenting organisms.^(2, 3) However very few organisms growing in the medium metabolize acetamide by the process of deamination (acrylamidase activity).^(4, 5) This unique ability is useful in identification of various non-fermenting gramnegative organisms.^(6, 7, 8) The ability to deaminate acetamide (acylaminidase activity) has been found to be most actively accomplished by *P. aeruginosa*, *Comamonas acidovorans*, *Achromobacter xylosoxidans* subsp. *xylosoxidans* (*Alcaligenes xylosoxidans*) and *Alcaligenes faecalis* (*odorans*).⁽⁹⁾ Acetamide deamination leads to the liberation of ammonia, which thereby increases the pH of the medium, leading to a subsequent colour change of the phenol red indicator from yellow orange to purplish red. Some strains require upto seven days to exhibit a positive reaction as they deaminate acrylamide slowly. However, only about 40% of apyocyanogenic strains of *Pseudomonas aeruginosa* exhibit a positive reaction. Therefore, it is advisable not to rely on this test as the only criterion for identification.

Principles of the Procedure

Acetamide Agar contains inorganic salts and acetamide a sole carbon and nitrogen source. Sodium chloride maintains the osmotic equilibrium. Phenol red is the pH indicator.

Formula / Liter

Ingredients	Gms / Liter
Part A	
Acetamide	10.00
Part B	
Sodium chloride	5.00
Dipotassium hydrogen phosphate	1.39
Potassium dihydrogen phosphate	0.73
Phenol red	0.012
Magnesium sulphate	0.50
Agar	15.00
Final pH: 7.0 ± 0.2 at 25°C	
Formula may be adjusted and/or supplemented as required to meet performance specifications	

Precautions

1. For Laboratory Use only.
2. IRRITANT. Irritating to eyes, respiratory system, and skin.

Directions

1. Suspend 22.63 grams of part B in 1000 ml distilled water.
2. Add 10.0 grams of Part A. Heat to boiling to dissolve the medium completely.
3. Dispense in tubes or as desired.
4. Autoclave at 121°C, 15 psi pressure, for 15 minutes / validated cycle.
5. Cool the tubes in a slanted position

Quality Control Specifications





PRODUCT SPECIFICATION SHEET

Dehydrated Appearance	Part A : Colourless deliquescent crystals Part B : Light yellow to brick red homogeneous free flowing powder
Prepared Medium	Orange coloured clear to slightly opalescent gel forms in tubes as slants
Reaction of (Mixture of 1% w/v Part A and 2.263% Part B) Solution	pH : 7.0 ± 0.2 at $25^\circ C$
Gel Strength	Firm, comparable with 1.5% Agar gel

Expected Cultural Response : Cultural characteristics observed after an incubation at $35-37^\circ C$ for 4-7 days.

Sr. No.	Organisms	Results to be achieved		
		Inoculum (CFU)	Growth	Deamination
1.	<i>Stenotrophomonas maltophilia</i> ATCC 13637	50-100	good- luxuriant	negative reaction, no purplish red colour within 7 days
2.	<i>Pseudomonas aeruginosa</i> ATCC 27853	50-100	good- luxuriant	positive reaction, purplish red colour within 7 days

The organisms listed are the minimum that should be used for quality control testing.

Test Procedure

1. Inoculate the Acetamide Agar slant with a loopful of culture emulsified in Trypticase Soy Broth.
2. Incubate inoculated slant at $35 \pm 2^\circ C$ and observe daily for 4 days and again at 7 days before discarding as negative.
3. Refer to appropriate references for standard test procedures.

Results

1. Deamination of the acetamide is indicated by a pronounced purplish-red color of the medium.
2. Complete identification requires determination of the Gram reaction, cellular morphology, biochemical reactions, etc.

Storage

Store the sealed bottle containing the dehydrated medium at $10 - 30^\circ C$. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light.

Expiration

Refer to the expiration date stamped on the container. The dehydrated medium should be discarded if not free flowing, or if the appearance has changed from the original color. Expiry applies to medium in its intact container when stored as directed.

Limitations of the Procedure

1. Some strains deaminate acetamide slowly and may require as long as 7 days to yield a positive test result.
2. Only about 37% of apyocyanogenic strains of *P. aeruginosa* will produce a positive reaction. Therefore, this test should not be relied upon as a sole criterion for identification.
3. For identification, organisms must be in pure culture. Morphological, biochemical and/or serological tests should be performed for final identification.
4. Consult appropriate texts for detailed information and recommended procedures.

Packaging

Product Name : Acetamide Agar (Twin Pack)

Product Code : DM444

Available Pack sizes : 500gm

References

1. Eaton A. D., Clesceri L. S., and Greenberg A. W., (Eds.), 1995, Standard Methods for the Examination of Water and Wastewater, 21st Ed., APHA, Washington, D.C.



PRODUCT SPECIFICATION SHEET



-
- 2. Gilardi, 1974, Antonie Van Leeuwenhoek, J. Microbiology Serol., 39:229.
 - 3. Stainer, Palleroni and Doudoroff, 1966, J. Gen Microbiol., 43:159.
 - 4. Pickett M. J. and Pedersen M.M., 1970, Can. J. Microbiol., 16:351.
 - 5. Pickett M. J. and Pedersen M.M., 1970, Can. J. Microbiol., 16:401.
 - 6. Hedberg, 1969, Appl. Microbiol., 17: 481
 - 7. Smith and Dayton, 1972, Appl. Microbiol., 24: 143
 - 8. Buhlmann, Vischer and Bruhin, 1961, J. Bacteriol., 82:787
 - 9. Oberhofer and Rowen, 1974, Appl. Microbiol., 28:720.

Further Information

For further information please contact your local MICROMASTER Representative.



MICROMASTER LABORATORIES PRIVATE LIMITED

DM444PSS,QAD/FR/024,Rev.00

Unit 38/39, Kalpataru Industrial Estate,
Off G.B. Road, Near 'R-Mall', Thane (W) - 400607. M.S. INDIA.
Ph: +91-9320126789/9833630009/9819991103
Email: sales@micromasterlab.com

Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

