

# PRODUCT SPECIFICATION SHEET



## GBS Selective Supplement (MS015)

### Intended Use

An Antibiotic Supplement, recommended for the selective isolation and presumptive identification of group B *Streptococci*.

### Composition

Each vial contains the following

*Ingredients	Concentration
Metronidazole	5.00mg
Gentamicin	1.00mg

### Directions

1. One vial, sufficient for 450ml medium.
2. Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water.
3. Mix well and add it to 475 ml sterile, molten, cooled (45-50°C) GBS Medium Base (DM115) along with 25 ml inactivated Horse Serum MS084.
4. Mix well.

### Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

**Note:** \* Not for medicinal use

### Further Information

For further information please contact your local MICROMASTER Representative.



**MICROMASTER LABORATORIES PRIVATE LIMITED**

**MS015PSS,Rev.00,Ver.00/01.02.2016**

Unit 38/39, Kalpataru Industrial Estate,  
Off G.B. Road, Near 'R-Mall', Thane (W) - 400607. M.S. INDIA.  
Ph: +91-22-25895505, 4760, 4681. Cell: 9320126789.  
Email: [micromaster@micromasterlab.com](mailto:micromaster@micromasterlab.com)

### Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.