



# PRODUCT SPECIFICATION SHEET

## CLEANROOM ENVIRONMENTAL MONITORING LOCATION MARKERS (CRT007G) Gamma Irradiated/CRT007(Non sterile)

### Intended Use

Steri-Master cleanroom environmental monitoring location markers used to easy identification of settle plate and air sampling locations in cleanroom area.

### Product Summary and Explanation

Steri-Master cleanroom environmental monitoring location markers are ideal for applications in and around the cleanroom. It is packed in a way so that it won't introduce any particles. Location markers is extremely conformable. The synthetic rubber adhesive ensures a 'clean peel' from most cleanroom work surfaces. It doesn't remains the residue. They provide effective moisture and vapor barrier. These location markers are chemically resistant and have a long shelf life. It is adhere to different surfaces

### Product Description

Bonds well with many metal, glass, cleanroom epoxy and plastic surfaces  
Ready to use

Backing	Vinyl
Adhesive	Rubber

### Precautions

1. Clean the area before applying
2. Area should be dust free and particle free
3. Surface should be smooth
4. Let it dry once cleaned
5. Peel off the sticker
6. Apply on your desired location
7. Apply pressure and make sure no air bubbles have been trapped and all peripheral area of the sticker should be stuck with surface.
8. Wait for at least 30min once applied. Do not make immediate cleaning.

### Quality Control Specifications

Appearance after Irradiation treatment:

Color	Change in color should not be observed.
Shape	Change in shape should not be observed.
Deterioration	Deterioration should not be observed.

The Vinyl tape comes with full certification of irradiation and lot traceability.

### Storage

Store at 15-30°C and 40 to 50% relative humidity.

### Expiration

To obtain best performance, use this product within 24 months from date of manufacture.





## PRODUCT SPECIFICATION SHEET

---

### Limitations of the Procedure

NA

### Packing Information

It is double layered packing containing 25/50 pieces of stickers. Each packet having Label & gamma irradiation indicator (For Sterile location markers)

### Packaging

**Product Name : Steri-Master Cleanroom Environmental Monitoring Location Markers**

**Product Code : CRT007**

**Available Pack sizes : As per customer requirement**

**Available Color Code : Black (BK), Red (RE), Yellow (YE)**

**Note: Packing and design of the location markers designed as per customers requirement.**

### Further Information

For further information please contact your local MICROMASTER Representative.



**MICROMASTER LABORATORIES PRIVATE LIMITED**

CRT007PSS,QAD/FR/024,Rev.00/01.01.2019

Unit 38/39, Kalpataru Industrial Estate,

Off G.B. Road, Near 'R-Mall', Thane (W) - 400607. M.S. INDIA.

Ph: +91-22-25895505, 4760, 4681. Cell: 9320126789.

Email: [micromaster@micromasterlab.com](mailto:micromaster@micromasterlab.com)

### Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.