



PRODUCT SPECIFICATION SHEET

Bromo Cresol Purple Azide Broth (DM387)

Intended Use

Bromo Cresol Purple Azide Broth (DM387) is recommended for the confirmatory presence of faecal *Streptococci* in water and waste water.

Product Summary and Explanation

Canned foods are most often prone to flat-sour spoilage due to contamination by either mesophilic or thermophilic aerobic spore-formers. In the 1930's, the National Canners Association specified use of Dextrose Tryptone Agar for isolating "flat sour" organisms from food products.⁽¹⁾ "Flat sour" spoilage of canned foods is caused by *Bacillus coagulans* (*Bacillus thermoacidurans*). Bacterial growth results in a 0.3-0.5 drop in pH, while ends of the can remain flat. *B.coagulans* is a soil microorganism, and found in canned tomato products and dairy products. Conditions favorable for multiplication of the bacterium can result in spoilage of the food product. Williams⁽²⁾ evolved Dextrose Tryptone Agar, a suitable medium for cultivation and enumeration of the thermophilic bacteria. It is also recommended for general cultural studies by Cameron⁽³⁾ and other associations.⁽⁴⁻⁸⁾ Dextrose Tryptone Broth, Modified is similar in composition to Dextrose Tryptone Agar, Modified, except agar. Dextrose Tryptone Broth, Modified is more nutritious and well buffered than Dextrose Tryptone Broth due to inclusion of yeast extract and dipotassium phosphate. This medium is useful for enumeration of mesophilic organisms, thermophiles in cereals and cereal products, dehydrated fruits and vegetables and spices.⁽⁹⁾

Principles of the Procedure

Dextrose Tryptone Broth Modified contains casein enzymic hydrolysate which provides carbon, nitrogen and other essential growth nutrients. Yeast extract supply essential vitamins. Dextrose serves as an energy source while bromocresol purple is a pH indicator. Dipotassium phosphate buffers the medium. Acid producing organisms produce yellow coloured medium. The tubes should be incubated at 55°C for 48 hours in a humid incubator.

Formula / Liter

Ingredients	Gms / Liter
Casein enzymic hydrolysate	10.00
Yeast extract	10.00
D-Glucose	5.00
Sodium chloride	5.00
Dipotassium hydrogen phosphate	2.70
Potassium dihydrogen phosphate	2.70
Sodium azide	0.50
Bromocresol purple	0.032
Final pH: 7.0 ± 0.2 at 25°C	
Formula may be adjusted and/or supplemented as required to meet performance specifications	

Precautions

1. For Laboratory Use only.
2. IRRITANT. Irritating to eyes, respiratory system, and skin.
3. Sodium azide has a tendency to form explosive metal azides with plumbing materials. It is advisable to use enough water to flush off the disposables.





PRODUCT SPECIFICATION SHEET

Directions

1. Suspend 35.93 grams of the medium in one liter of distilled water. Add 5 ml glycerol if desired.
2. Heat if necessary to dissolve the medium completely.
3. Dispense into test tubes.
4. Autoclave at 121°C, 15 psi pressure, for 15 minutes / validated cycle.

Quality Control Specifications

Dehydrated Appearance	Light yellow to beige homogeneous free flowing powder
Prepared Medium	Purple coloured, clear solution without any precipitate
Reaction of 3.6% Solution	pH : 7.0 ± 0.2 at 25°C
Gel Strength	Not Applicable

Expected Cultural Response: Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Sr. No.	Organisms	Results to be achieved		
		Inoculum (CFU)	Growth	Acid
1.	<i>Enterococcus faecalis</i> ATCC 29212	50-100	good-luxuriant	positive reaction, yellow colour
2.	<i>Escherichia coli</i> ATCC 25922	>=10 ³	inhibited	--
3.	<i>Pseudomonas aeruginosa</i> ATCC 27853	>=10 ³	inhibited	--
4.	<i>Streptococcus agalactiae</i> ATCC 13813	50-100	none-poor	negative reaction, no colour change
5.	<i>Streptococcus pyogenes</i> ATCC 19615	50-100	none-poor	negative reaction, no colour change

The organisms listed are the minimum that should be used for quality control testing.

Test Procedure

Refer appropriate references for standard test procedures.

Results

Refer appropriate references and procedures for interpretation of results.

Storage

Store the sealed bottle containing the dehydrated medium at 10 - 30°C. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light.

Expiration

Refer to the expiration date stamped on the container. The dehydrated medium should be discarded if not free flowing, or if the appearance has changed from the original color. Expiry applies to medium in its intact container when stored as directed.

Limitations of the Procedure





PRODUCT SPECIFICATION SHEET

1. For identification, organisms must be in pure culture. Morphological, biochemical and/or serological tests should be performed for final identification.
2. Consult appropriate texts for detailed information and recommended procedures.

Packaging

Product Name : Bromo Cresol Purple Azide Broth

Product Code : DM387

Available Pack sizes : 500gm

References

1. National Canners Association. 1933. Bacterial standards for sugar.
2. Williams O.B., 1936, Food Res., 1:217.
3. Cameron E.J., 1936, J.Assoc. Official Agr. Chem., 19:433.
4. Association of Official Analytical Chemists, 1978, Bacteriological Analytical Manual, 5th Edition, AOAC, Washington, D.C.
5. American Public Health Association, 1972, Standard Methods for the Examination of Dairy Products, 13th Ed. APHA, Washington, D.C.
6. National Canners Association, 1968, Laboratory Manual for Food Caners and Processors, Vol. I.
7. American Public Health Association, 1976, Compendium of Methods for the Microbiological Examination of Foods, APHA, Washington, D.C.
8. National Canners Association, 1954, A Laboratory Manual for the Canning Industry, 1st Edition, National Canners Associations, Washington.
9. Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.

Further Information

For further information please contact your local MICROMASTER Representative.



MICROMASTER LABORATORIES PRIVATE LIMITED




DM387PSS,QAD/FR/024,Rev.00/01.01.2018

Unit 38/39, Kalpataru Industrial Estate,
Near Runwal Estate, Behind 'R-Mall' ,Ghodbunder Raod,
Thane (W) - 400607. M.S. INDIA.

Ph: +91-22-25895505, 4760, Cell: 9320126789.

Email: micromaster@micromasterlab.com

sales@micromasterlab.com

Prepared By	Checked By	Approved By
 01.01.2018	 01.01.2018	 01.01.2018
Microbiologist	Head Quality Control	Head Quality Assurance





PRODUCT SPECIFICATION SHEET

Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

