



Sterimaster

CLEANROOM ACCESSORIES



Common features:

Processing and Packing in cleanroom environment in adherence ISO Class 4 and EUGMP Grade A

Gamma Irradiated triple-bagged packaging that reduces particles

Supportive Real-time validation studies

Ready to use

Complete traceability

CLEANROOM STATIONERY

1. Cleanroom pen and Marker:

Non-smearing, low odor ink

2. Cleanroom Paper:

Compatible with regular ink or photocopying

Advantages

Coated with specialized polymer


Autoclavable

Intangible benefit

Durable,

Does not generate particles

Cost-effective and saves time



CLEANROOM TAPES

Advantages

Specially designed core

Synthetic rubber adhesive

Synthetic material

Autoclavable

No requirement of Aluminum foil covers

Firm Grip

Intangible benefit

Does not generate particles

No yeast and mould growth

No residue

Adhere to Different surfaces

Written on with pen and marker

Cost-effective and saves time

Increases aesthetic and visibility of packed plate

Decreases risk of plates falling or opening

Accidental hazards and false positive

results avoided



FLOOR MARKING TAPES

Advantages

Specially designed core

Synthetic rubber adhesive

No specialized skill required to demark the area

Resistant to peeling and disinfectants

Intangible benefit

Does not generate particles in cleanroom

No yeast and mould growth


No residue

Adhere to Different surfaces

Easy to use

Cost effective and saves time

Repeated painting costs are avoided



We also have:

STERILE SWAB

STERILE FLEXILOOP

AUTOCLAVABLE BAGS

CLEANROOM BOTTLES

