

# CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY



Product Name: Brain Heart Infusion Broth	
Product Code: DM811	Expiry Date: 06/02/2021
Lot Number: 00630116	Release Date: 06/02/2016
Analytical Report Number: QAD/AR-15/02-244	Page 1 of 2

DEHYDRATED MEDIUM			
TEST CONDUCTED	SPECIFICATION	OBSERVATIONS	REMARKS
Color	Cream to yellow colored powder	Cream to yellow colored powder	Complies
Consistency	Homogenous, free flowing powder	Homogenous, free flowing powder	Complies
PREPARED MEDIUM			
TEST CONDUCTED	SPECIFICATION	OBSERVATIONS	REMARKS
Color of medium	Light to medium amber colored solution	Light to medium amber colored solution	Complies
Clarity of medium (After sterilization)	Clear solution without any precipitate.	Clear solution without any precipitate.	Complies
pH (After sterilization)	7.20 to 7.60 at 25°C	7.50	Complies
Gel strength	Not Applicable	Not Applicable	Complies

Rev. 00


10/12/2014

QAD/FR/01

**Micromaster Laboratories Pvt. Ltd.**  
Manufacture & Exporter

Culturing your needs in Microbiology

38/39, Kalpataru Industrial Estate, Opp Lawkim, Off. Ghodbunder Rd.,  
Near Runwal Estate, Thane (W) 400 607, Maharashtra, INDIA.  
Ph: +91-22-2589 5505, 4760, 4681, 9320126789  
Email: micromaster@micromasterlab.com



# CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY



Product Name: Brain Heart Infusion Broth	
Product Code: DM811	Expiry Date: 06/02/2021
Lot Number: 00630116	Release Date: 06/02/2016
Analytical Report Number: QAD/AR-15/02-244	Page 2 of 2

## CULTURAL PARAMETERS:

Organisms	Incubation time/ temp	No. of organisms inoculated	Growth	Remarks
<i>Enterococcus faecalis</i> ATCC 29212	35-37°C / 24-48 hrs.	84	Good-Luxuriant	Complies
<i>Neisseria meningitidis</i> ATCC 13090	35-37°C / 24-48 hrs.	65	Good-Luxuriant	Complies
<i>Streptococcus pneumoniae</i> ATCC 6303	35-37°C / 24-48 hrs.	68	Good-Luxuriant	Complies
<i>Streptococcus pyogenes</i> ATCC 19615	35-37°C / 24-48 hrs.	64	Good-Luxuriant	Complies
<i>Candida albicans</i> ATCC 10231	35-37°C / 24-48 hrs.	72	Good-Luxuriant	Complies
<i>Staphylococcus aureus</i> ATCC 25923	35-37°C / 24-48 hrs.	93	Good-Luxuriant	Complies

## Final Results: COMPLIES

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

TESTED BY (Vijayalaxmi Pujari) MICROBIOLOGIST	CHECKED BY (Swati More) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report

Rev. 00

10/12/2014

QAD/FR/01

Micromaster Laboratories Pvt. Ltd.  
Manufacture & Exporter

Culturing your needs in Microbiology

38/39, Kalpataru Industrial Estate, Opp Lawkim, Off. Ghodbunder Rd.,  
Near Rurwal Estate, Thane (W) 400 607, Maharashtra, INDIA.  
Ph: +91-22-2589 5505, 4760, 4661, 9320126789  
Email: micromaster@micromasterlab.com

