

CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY



Product Name: Sterile Specimen Collector	
Product Code: PW001-30ml	Expiry Date: 10/12/2020
Lot Number: 01361115	Release Date: 10/12/2015
Analytical Report Number: QAD/AR-15/12-302	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Sterilized specimen Collector with Individually sealed, pack of 100	Sterilized specimen Collector with Individually sealed, pack of 100.	Complies
Sterilization	The above batch of specimen Collector was sterilized by Microtrol Sterilization Services Pvt. Ltd. by treatment with Ethylene oxide to the following parameters: Concentration of Ethylene oxide : 750 gm/M ³ Exposure time : 6hours Humidity : Min. 50% RH Temperature : 50°C	Concentration of Ethylene oxide: 750 gm/M ³	Complies
Sterility Testing	As per USP33/NF28.	Sterility passes	Complies

Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

TESTED BY (Sneha Daragshetti) MICROBIOLOGIST	CHECKED BY (Swati More) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report

Rev. 00

10/12/2014

QAD/FR/01

Micromaster Laboratories Pvt. Ltd.
Manufacture & Exporter

Culturing your needs in Microbiology...

38/39, Kalpataru Industrial Estate, Opp Lawkim, Off. Ghodbunder Rd.,
Near Runwal Estate, Thane (W) 400 607, Maharashtra, INDIA.
Ph: +91-22-2589 5505, 4760, 4681, 9320126789
Email: micromaster@micromasterlab.com

