

PRODUCT SPECIFICATION SHEET



Potassium Tellurite Supplement-3.5% (MS025)(1ml per vial)

Intended Use

A Selective Supplement, recommended for selective isolation of *Staphylococci* and *Corynebacteria*.

Composition

Each vial contains the following:

*Ingredients	Concentration
Potassium Tellurite	0.35 gm
Distilled water	1.00 ml

Directions

1. One vial, sufficient for 1000 ml sterile molten medium.
2. Warm up the refrigerated 3.5% Sterile Potassium Tellurite Solution (MS025) to 45-50° C. Add aseptically 9 ml sterile distilled water, mix well and add 3 ml in 1000 ml Baird parker Agar Base (DM034) or 10 ml in 1000ml Hoyle Medium Base (DM120) or 2 ml in 1000 ml Cholera Medium Base (DM314).
3. Mix gently and pour into sterile petri plates.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Note: * Not for medicinal use

Further Information

For further information please contact your local MICROMASTER Representative.



MICROMASTER LABORATORIES PRIVATE LIMITED

MS025PI, Rev.00, Ver.01/02.12.2016

Unit 38/39, Kalpataru Industrial Estate,
Off G.B. Road, Near 'R-Mall', Thane (W) - 400607. M.S. INDIA.
Ph: +91-22-25895505, 4760, 4681. Cell: 9320126789.
Email: micromaster@micromasterlab.com

Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.