






CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: Master anaerobic system PG Jar	
Product Code: AAN004	Expiry Date: Not applicable
Lot Number: 1253	Release Date: 24/08/2018
Analytical Report Number: QAD/AR-18/08-803	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATION	OBSERVATIONS	REMARKS
Description	Polycarbonate Jar	Polycarbonate Jar	Complies
Range	30 psi (-1Kg/cm ²) to +30 psi (2.1Kg/cm ²)	30 psi (-1Kg/cm ²) to +30 psi (2.1Kg/cm ²)	Complies
Capacity	3.5lit	3.5lit	Complies
Use	Used to create anaerobic condition for the growth of anaerobic microorganism		

Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
TESTED BY (Madhura Kadam) MICROBIOLOGIST	CHECKED BY (Aditi Modak) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report