






## CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

<b>Product Name: Master anaerobic system PG Jar</b>	
<b>Product Code: AAN004</b>	<b>Expiry Date: Not applicable</b>
<b>Lot Number: 1254</b>	<b>Release Date: 24/08/2018</b>
<b>Analytical Report Number: QAD/AR-18/08-803</b>	<b>Page No: 1 of 1</b>

TEST CONDUCTED	SPECIFICATION	OBSERVATIONS	REMARKS
Description	Polycarbonate Jar	Polycarbonate Jar	Complies
Range	30 psi (-1Kg/cm <sup>2</sup> ) to +30 psi ( 2.1Kg/cm <sup>2</sup> )	30 psi (-1Kg/cm <sup>2</sup> ) to +30 psi ( 2.1Kg/cm <sup>2</sup> )	Complies
Capacity	3.5lit	3.5lit	Complies
Use	Used to create anaerobic condition for the growth of anaerobic microorganism		

### Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
<b>TESTED BY (Madhura Kadam) MICROBIOLOGIST</b>	<b>CHECKED BY (Aditi Modak) HEAD MICROBIOLOGY</b>	<b>APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE</b>

End Of The Report