

# CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY






<b>Product Name: Ehrlich Reagent</b>	
<b>Product Code: AR020</b>	<b>Expiry Date : 18/04/2025</b>
<b>Lot Number: 02500417</b>	<b>Release Date: 18/04/2017</b>
<b>Analytical Report Number: QAD/AR-17/04-439</b>	<b>Page No : Page 1 of 1</b>

TESTS CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Light brown colored	Light brown colored	Complies.
Use	Used for to detect urobilinogen in fresh, cool urine		

## Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
<b>TESTED BY (Daniya Kazi) MICROBIOLOGIST</b>	<b>CHECKED BY (Vijayalakshmi Poojari) HEAD MICROBIOLOGY</b>	<b>APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE</b>

End Of The Report

Rev. 00

10/12/2014

QAD/FR/01

**Micromaster Laboratories Pvt. Ltd.**  
Manufacture & Exporter

Culturing your needs in Microbiology...

38/39, Kalpataru Industrial Estate, Opp Lawkim, Off. Ghodbunder Rd.,  
Near Runwal Estate, Thane (W) 400 607, Maharashtra, INDIA.  
Ph: +91-22-2589 5505, 4760, 4681, 9320126789  
Email: micromaster@micromasterlab.com

