






CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: Cleanroom Spray Bottle	
Product Code: CRA001	Expiry Date: Not Applicable
Lot Number: 03260618	Release Date: 04/06/2018
Analytical Report Number: QAD/AR-18/06-520	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Semi transparent , Polypropylene, screw cap bottle with trigger , Individually packed		
Steam Sterilization (121°C at 15 psi)	Bottle should withstand steam Sterilization (121°C at 15 psi)-After sterilization no any change should observed.	Bottle withstands steam Sterilization (121°C at 15 psi)-After sterilization no any change was observed	Complies

Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
TESTED BY (Madhura Kadam) MICROBIOLOGIST	CHECKED BY (Aditi Modak) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report