



CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY




| | |
|---|------------------------------------|
| Product Name: Cleanroom Spray Bottle | |
| Product Code: CRA001 | Expiry Date: Not applicable |
| Lot Number: G09101718 | Release Date: 22/07/2018 |
| Analytical Report Number: QAD/AR-18/07-670 | Page No: 1 of 1 |

| TEST CONDUCTED | SPECIFICATIONS | OBSERVATIONS | REMARKS |
|--|---|---|----------|
| Appearance | Semi transparent, Polypropylene, screw cap bottle with trigger. Packing-Eight Spray Head and one bottle | | |
| Steam Sterilization (121°C at 15 psi) | Bottle and spray Head should withstand steam Sterilization (121°C at 15 psi)-After sterilization no any change should observed. | Bottle and spray Head withstands steam Sterilization (121°C at 15 psi)-After sterilization no any change was observed | Complies |

Remark: Spray Head- Autoclavable- Upto 5 Cycles
Bottle- Autoclavable- Upto 40 Cycles

Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

| | | |
|---|---|---|
|  |  |  |
| TESTED BY Madhura Kadam (MICROBIOLOGIST) | CHECKED BY Aditi Modak (HEAD MICROBIOLOGY) | APPROVED BY Rupali Sabale (HEAD QUALITY ASSURANCE) |

End Of The Report

CRA001COA,QAD/FR/031,Rev.00/21.07.2018

Micromaster laboratories Pvt.Ltd

Manufacturer and Exporter

Unit No.38/39, Kalpataru Industrial Estate .Near Runwal Estate,Behind R Mall,
Ghodbunder Road ,Thane (W)-400607,Maharashtra,India

Tel:- +91-22-2589 5505sales@micromasterlab.com