

CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY



Product Name: Master-Wash bottle	
Product Code: CRA002	Expiry Date: NA
Lot Number: 03740617	Release Date: 09/06/2017
Analytical Report Number: QAD/AR-17/06-649	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Semi transparent , Polypropylene, screw cap bottle with nozzle & cap , Individually packed	Semi transparent , Polypropylene, screw cap bottle with nozzle & cap , Individually packed bottle was observed	Complies
Steam Sterilization (121°C at 15 psi)	Bottle should withstand steam Sterilization (121°C at 15 psi)	Bottle withstands steam Sterilization (121°C at 15 psi)	Complies

Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

TESTED BY (Divya Sasidharan) MICROBIOLOGIST	CHECKED BY (Vijayalakshmi Poojari) HEAD MICROBIOLOGY	APPROVED BY (Tejal Shinde) HEAD QUALITY ASSURANCE

End Of The Report

Rev. 00

10/12/2014

QAD/FR/01

Micromaster Laboratories Pvt. Ltd

Manufacturer & Exporter

Unit No.38/39, Kalpataru Industrial Estate, Near Runwal Estate, Behind R Mall,
Ghodbunder Road, Thane (W) - 400607, Maharashtra, India

Tel : +91-22-2589 5505 / 2589 4681 / 2589 4760 sales@micromasterlab.com

