






CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: CleanRoom Reagent Bottle	
Product Code: CRA004	Expiry Date: NA
Lot Number: 03250618	Release Date: 14/06/2018
Analytical Report Number: QAD/AR-18/06-557	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Semi transparent , Polypropylene, screw cap bottle, Individually packed	Semi transparent , Polypropylene, screw cap bottle , Individually packed bottle was observed	Complies
Steam Sterilization (121°C at 15 psi)	The above batch of clean room reagent bottle was sterilized by ISOMED (BriT) by treatment with Gamma irradiation to following parameters: Dose of irradiation 10-25 kgy.	Dose of irradiation 16.9 kgy	Complies

Final Results: COMPLIES

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
TESTED BY (Madhura Kadam) MICROBIOLOGIST	CHECKED BY (Aditi Modak) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report