



CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: Sterile Transport Swab(Swab w/Tube) 8"	
Product Code: CRA007-100PC	Expiry Date: 06/11/2024
Lot Number: 05130619	Release Date: 13/07/2019
Analytical Report Number: QAD/AR-19/07-0731	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Apearance	Transparent, clear & unbreakable, polypropylene tube with specially designed cap and cotton bud with polypropylene stick. Individually packed, Size: 150mm x 12mm diameter. Pack of 100 nos.	Transparent, clear & unbreakable, polypropylene tube with specially designed cap and cotton bud with polypropylene stick. Individually packed, Size: 150mm x 12mm diameter. Pack of 100 nos.	Complies
Sterilization	The above batch of swabs was sterilized by Microtrol Sterilization Services Pvt. Ltd. by treatment with Ethylene oxide to the following parameters: Concentration of Ethylene oxide : NLT 750 gm/M3 Exposure time : 6 hours Humidity : Min. 50% RH Temperature : 50 Deg C	Concentration of Ethylene oxide: 750 gm/M3	Complies
Sterility Testing	Test result should comply with sterility test (71) of USP 37 method.	Sterility passes	Complies

Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
TESTED BY (Pooja Sutar) MICROBIOLOGIST	CHECKED BY (Madhavi Sawant) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report

CRA007COA,QAD/FR/031,Rev.00/01.01.2018

Micromaster laboratories Pvt.Ltd
Manufacturer and Exporter
Unit No.38/39, Kalpataru Industrial Estate .Near Runwal Estate,Behind R Mall,
Ghodbunder Road,Thane (W)-400607,Maharashtra,India
[Tel:- +91-22-2589 5505sales@micromasterlab.com](mailto:+91-22-2589 5505sales@micromasterlab.com)