






## CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

<b>Product Name: Autoclavable Petri plates,Transparent-90mm</b>	
<b>Product Code: CRA016-10PC</b>	<b>Expiry Date: 19/07/2023</b>
<b>Lot Number: G0900718</b>	<b>Release Date: 14/07/2018</b>
<b>Analytical Report Number: QAD/AR-18/07-627</b>	<b>Page No: 1 of 1</b>

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Polycarbonate, clear, transparent & Unbreakable.	Polycarbonate, clear, transparent & Unbreakable.	Complies
Size	90mm diameter x 15mm	90mm diameter x 15mm	Complies
Steam Sterilization (121°C at 15 psi)	Bottle should withstand steam Sterilization (121°C at 15 psi)	Bottle withstands steam Sterilization(121°C at 15 psi)	Complies
Use	Used for Microbiological testing.		

### **Final Results: COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
<b>TESTED BY (Madhura Kadam) MICROBIOLOGIST</b>	<b>CHECKED BY (Aditi Modak) HEAD MICROBIOLOGY</b>	<b>APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE</b>

End Of The Report