






## CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: Sterile Transparent Self Locking Poly Bags 6" x 9"	
Product Code: CRA027-500PC	Expiry Date: 22/07/2023
Lot Number: 04060718	Release Date: 30/07/2018
Analytical Report Number: QAD/AR-18/07-645	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Transparent, Self Locking Size: 6" x 9"	Transparent, Self Locking Size: 6" x 9"	Complies
Sterilization	The above batch of clean room tape was sterilized by AVI PROCESSOR by treatment with Gamma irradiation to the following parameters:  Dose of irradiation 10-25 kgy		
After Gamma Irradiation	1.No color change	1.No color change	Complies
	2.No any deterioration	2.No any deterioration	Complies

### Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
<b>TESTED BY (Madhura Kadam) MICROBIOLOGIST</b>	<b>CHECKED BY (Aditi Modak) HEAD MICROBIOLOGY</b>	<b>APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE</b>

End Of The Report