






CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: Sterile Transparent Self Locking Poly Bags 10"x12"	
Product Code: CRA028	Expiry Date: 27/07/2023
Lot Number: 04500718	Release Date: 03/08/2018
Analytical Report Number: QAD/AR-18/08-727	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Transparent,Self Locking Size: 10"x12"	Transparent,Self Locking Size: 10"x12"	Complies
Sterilization	The above batch of clean room tape was sterilized by AVI PROCESSOR by treatment with Gamma irradiation to the following parameters: Dose of irradiation 10-25 kgy		
After Gamma Irradiation	1.No color change	1.No color change	Complies
	2.No any deterioration	2.No any deterioration	Complies

Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
TESTED BY (Madhura Kadam) MICROBIOLOGIST	CHECKED BY (Aditi Modak) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report

CRA028COA,QAD/FR/031,Rev.00/01.01.2018

Micromaster laboratories Pvt.Ltd

Manufacturer and Exporter

Unit No.38/39, Kalpataru Industrial Estate .Near Runwal Estate,Behind R Mall,
Ghodbunder Road ,Thane (W)-400607,Maharashtra,India

Tel:- +91-22-2589 5505sales@micromasterlab.com