






## CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

<b>Product Name: Sterile Transparent Self Locking Bags 4x6"</b>	
<b>Product Code: CRA033</b>	<b>Expiry Date: 27/07/2023</b>
<b>Lot Number: 04510718</b>	<b>Release Date: 03/08/2018</b>
<b>Analytical Report Number: QAD/AR-18/08-728</b>	<b>Page No: 1 of 1</b>

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Transparent, Self Locking Size: 4"x6"	Transparent, Self Locking Size: 4"x6"	Complies
Sterilization	The above batch of clean room tape was sterilized by AVI PROCESSOR by treatment with Gamma irradiation to the following parameters:  Dose of irradiation 10-25 kgy		
After Gamma Irradiation	1.No color change	1.No color change	Complies
	2.No any deterioration	2.No any deterioration	Complies

### Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
<b>TESTED BY (Madhura Kadam) MICROBIOLOGIST</b>	<b>CHECKED BY (Aditi Modak) HEAD MICROBIOLOGY</b>	<b>APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE</b>

End Of The Report

CRA033COA,QAD/FR/031,Rev.00/01.01.2018

**Micromaster laboratories Pvt.Ltd**

Manufacturer and Exporter

Unit No.38/39, Kalpataru Industrial Estate .Near Runwal Estate,Behind R Mall,

Ghodbunder Road ,Thane (W)-400607,Maharashtra,India

Tel:- +91-22-2589 5505sales@micromasterlab.com