






# CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

<b>Product Name: Steri-Master Marker-Black</b>	
<b>Product Code: CRS002-BK-G</b>	<b>Expiry Date: 10/06/2019</b>
<b>Lot Number: 03600617</b>	<b>Release Date: 19/06/2017</b>
<b>Analytical Report Number: QAD/AR-17/06-629</b>	<b>Page No: 1 of 1</b>

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Individually packed, Black permanent, quick drying, alcohol resistant ink	Individually packed, Black permanent, quick drying, alcohol resistant ink	Complies
Sterilization	The above batch of Clean Room Marker was sterilized by ISOMED (BRIT) by irradiation to the following parameters:  Dose of irradiation 10-25 kgy.	Dose of irradiation 15Kgy	Complies
Sterility Testing	As per USP36/NF31.	Sterility passes	Complies

## Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
<b>TESTED BY (Tejal Shinde) MICROBIOLOGIST</b>	<b>CHECKED BY (Vijayalakshmi Poojari) HEAD MICROBIOLOGY</b>	<b>APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE</b>

End Of The Report

Rev. 00

10/12/2014

QAD/FR/01

**Micromaster Laboratories Pvt. Ltd**

Manufacturer & Exporter

Unit No.38/39, Kalpataru Industrial Estate, Near Runwal Estate, Behind R Mall,  
Ghodbunder Road, Thane (W) - 400607, Maharashtra, India

Tel : +91-22-2589 5505 / 2589 4681 / 2589 4760 sales@micromasterlab.com

