






## CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

<b>Product Name: Cleanroom Vinyl Tape White</b>	
<b>Product Code: CRT001WHG</b>	<b>Expiry Date: 05/03/2020</b>
<b>Lot Number: 01000218</b>	<b>Release Date: 08/03/2018</b>
<b>Analytical Report Number: QAD/AR-18/02-172</b>	<b>Page No : Page 1 of 1</b>

TESTS CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Color	White	White	Complies
Thickness	6.0 ± 0.7 mil/ 0.15 ± 0.02 mm	6.2 mil/ 0.157 mm	Complies
Sterilization	The above batch of clean room tape was sterilized by ISOMED (BRIT) by treatment with Gamma irradiation to the following parameters: Dose of irradiation 10-25 kgy		
After Gamma Irradiation	1.No color change	1.No color change	Complies
	2.No any deterioration	2.No any deterioration	Complies

### Final Results: COMPLIES

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
<b>TESTED BY (Deepti Nair) MICROBIOLOGIST</b>	<b>CHECKED BY (Aditi Modak) HEAD MICROBIOLOGY</b>	<b>APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE</b>

End of the Report