



CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: Diluting Fluid K	
Product Code: DM1612	Expiry Date: 14/04/2022
Lot Number: 02220417	Release Date: 14/04/2017
Analytical Report Number: QAD/AR-17/04-412	Page No: 1 of 2

DEHYDRATED MEDIUM			
TEST CONDUCTED	SPECIFICATION	OBSERVATIONS	REMARKS
Color	Cream to yellow colored powder	Cream to yellow colored powder	Complies
Consistency	Homogenous, free flowing powder	Homogenous, free flowing powder	Complies
PREPARED MEDIUM			
TEST CONDUCTED	SPECIFICATION	OBSERVATIONS	REMARKS
Color of medium	Light yellow colored medium	Light yellow colored medium	Complies
Clarity of medium (after sterilization)	Clear solution	Clear solution	Complies
pH (after sterilization)	6.70-7.10 at 25°C	6.77	Complies
Gel strength	Not Applicable	Not Applicable	Not Applicable

CULTURAL PARAMETERS: Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organisms	Incubation Time/ temp	No. of organisms inoculated	Growth	Remarks
<i>Candida albicans</i> ATCC 10231	35-37°C for 24-48 hours.	74	Good	Complies
<i>Escherichia coli</i> ATCC 25922	35-37°C for 24-48 hours.	80	Good	Complies
<i>Staphylococcus aureus</i> ATCC 25923	35-37°C for 24-48 hours.	81	Good	Complies

Rev. 00

10/12/2014

QAD/FR/01

Micromaster Laboratories Pvt. Ltd

Manufacturer & Exporter

Unit No.38/39, Kalpataru Industrial Estate, Near Runwal Estate, Behind R Mall,
Ghodbunder Road, Thane (W) - 400607, Maharashtra, India

Tel : +91-22-2589 5505 / 2589 4681 / 2589 4760 sales@micromasterlab.com





CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY




Product Name: Diluting Fluid K	
Product Code: DM1612	Expiry Date: 14/04/2022
Lot Number: 02220417	Release Date: 14/04/2017
Analytical Report Number: QAD/AR-17/04-412	Page No: 2 of 2

CULTURAL PARAMETERS: Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organisms	Incubation Time/ temp	No. of organisms inoculated	Growth	Remarks
<i>Escherichia coli</i> ATCC 8739	35-37°C for 24-48 hours	79	Good	Complies
<i>Staphylococcus aureus</i> ATCC 6538	35-37°C for 24-48 hours	77	Good	Complies

Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
TESTED BY (Daniya Kazi) MICROBIOLOGIST	CHECKED BY (Vijayalakshmi Poojari) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report

Rev. 00

10/12/2014

QAD/FR/01

Micromaster Laboratories Pvt. Ltd

Manufacturer & Exporter

Unit No.38/39, Kalpataru Industrial Estate, Near Runwal Estate, Behind R Mall,
Ghodbunder Road, Thane (W) - 400607, Maharashtra, India

Tel : +91-22-2589 5505 / 2589 4681 / 2589 4760 sales@micromasterlab.com

