

CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY



Product Name: Asparagine Proline Broth	
Product Code: DM318	Expiry Date: 02/02/2022
Lot Number: 00580117	Release Date: 02/02/2017
Analytical Report Number: QAD/AR-17/02-204	Page No: 1 of 1

DEHYDRATED MEDIUM			
TEST CONDUCTED	SPECIFICATION	OBSERVATIONS	REMARKS
Color	White to cream colored powder	White to cream colored powder	Complies
Consistency	Homogenous, free flowing powder	Homogenous, free flowing powder	Complies
PREPARED MEDIUM			
TEST CONDUCTED	SPECIFICATION	OBSERVATIONS	REMARKS
Color of medium	Colourless clear Solution	Colourless clear Solution	Complies
Clarity of medium (After sterilization)	Clear solution without any precipitate	Clear solution without any precipitate	Complies
pH (After sterilization)	7.10 to 7.50 at 25°C	7.28	Complies
Gel strength	Not Applicable	Not Applicable	Not Applicable

CULTURAL PARAMETERS:

Organisms	Incubation time/ temp	No. of organisms inoculated	Growth	Remarks
<i>Escherichia coli</i> ATCC25922	35-37°C / 24-48 hrs.	89	None to poor	Complies
<i>Pseudomonas aeruginosa</i> ATCC 27853	35-37°C / 24-48 hrs.	84	Luxuriant with greenish yellow pigment	Complies

Final Results: COMPLIES

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

TESTED BY (Divya Sasidharan) MICROBIOLOGIST	CHECKED BY (Tejal Shinde) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report

Rev. 2

01/10/2013

QAD/FR/01

Micromaster Laboratories Pvt. Ltd.
Manufacture & Exporter

Culturing your needs in Microbiology...

38/39, Kalpataru Industrial Estate, Opp Lawkim, Off. Ghodbunder Rd.,
Near Runwal Estate, Thane (W) 400 607, Maharashtra, INDIA.
Ph: +91-22-2589 5505, 4760, 4681, 9320126789
Email: micromaster@micromasterlab.com

