



PRODUCT SPECIFICATION SHEET

Tributyryn (10ml per vial) (MS046)

Intended Use

A Supplement, recommended for detection of lipolytic microorganisms.

Composition

Each vial contains the following

| *Ingredients | Concentration |
|--------------|---------------|
| Tributyryn | 10ml |

Directions

1. One vial, sufficient for 990ml medium.
2. Warm up the refrigerated solution to room temperature.
3. Add contents of one vial aseptically to 990 ml of Tributyrin Agar Base w/o Tributyrin (DM270) and then sterilize by autoclaving at 15 lbs pressure for 15 minutes.
4. Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Note: * Not for medicinal use

Further Information

For further information please contact your local MICROMASTER Representative.



MICROMASTER LABORATORIES PRIVATE LIMITED

MS046PSS,QAD/FR/024,Rev.00/01.01.2018

Unit 38/39, Kalpataru Industrial Estate,
Off G.B. Road, Near 'R-Mall', Thane (W) - 400607. M.S. INDIA.
Ph: +91-22-25895505, 4760, 4681. Cell: 9320126789.
Email: micromaster@micromasterlab.com

| Prepared By | Checked By | Approved By |
|---|---|---|
|  01.01.2018 |  01.01.2018 |  01.01.2018 |
| Microbiologist | Head Quality Control | Head Quality Assurance |





PRODUCT SPECIFICATION SHEET

Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

