

PRODUCT SPECIFICATION SHEET



Kimmig Selective Supplement (Twin Pack) (MS081)

Intended Use

An Antibiotic Supplement, recommended for selective isolation and cultivation of fungi.

Composition

Each vial contains the following

*Ingredients	Concentration
Part A	
Colistin sulphate	40.00mg
Part B	
Novobiocin	50.00mg

Directions

1. One vial, sufficient for 500ml medium.
2. Rehydrate the contents of one vial each of Part A and Part B, separately with 2.5 ml sterile distilled water and mix well.
3. On reconstitution aseptically add to 500 ml sterile, molten, cooled (45-50°C) Kimmig Fungi Agar Base (DM1390).
4. Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

Note: * Not for medicinal use

Further Information

For further information please contact your local MICROMASTER Representative.



MICROMASTER LABORATORIES PRIVATE LIMITED

MS081PSS,Rev.00,Ver.00/01.02.2016

Unit 38/39, Kalpataru Industrial Estate,
Off G.B. Road, Near 'R-Mall', Thane (W) - 400607. M.S. INDIA.
Ph: +91-22-25895505, 4760, 4681. Cell: 9320126789.
Email: micromaster@micromasterlab.com

Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.