



PRODUCT SPECIFICATION SHEET

L. mono Selective Supplement-I (MS157)

Intended Use

A Selective Supplement, recommended for isolation of *Listeria* species.

Composition

Each vial contains the following:

*Ingredients	Concentration
Polymyxin B sulphate	38350IU

Directions

1. One vial, sufficient 500ml medium.
2. Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water.
3. Mix well and aseptically add it to 460 ml of sterile, molten, cooled (45-50°C) L. mono Differential Agar Base (DM1543) along with sterile contents of one vial of L. mono Enrichment Supplement I (MS067) and sterile rehydrated contents of one vial of L. mono Selective Supplement II (MS158) or add in 470 ml of sterile, molten, cooled (45-50°C) L. mono Confirmatory Agar Base (DM1538) along with sterile contents of one vial of L. mono Enrichment Supplement II (MS159) and rehydrated contents of one vial of L. mono Selective Supplement II (MS158).
4. Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2- 8°C. Use before expiry date on the label.

Note: * Not for medicinal use

Further Information

For further information please contact your local MICROMASTER Representative.



MICROMASTER LABORATORIES PRIVATE LIMITED




MS157PSS, QAD/FR/024, Rev.00/01.01.2018

Unit 38/39, Kalpataru Industrial Estate,
Off G.B. Road, Near 'R-Mall', Thane (W) - 400607. M.S. INDIA.
Ph: +91-22-25895505, 4760, 4681. Cell: 9320126789.
Email: micromaster@micromasterlab.com





PRODUCT SPECIFICATION SHEET

Prepared By	Checked By	Approved By
 01.01.2018	 01.01.2018	 01.01.2018
Microbiologist	Head Quality Control	Head Quality Assurance

Disclaimer :

All Products conform exclusively to the information contained in this and other related Micromaster Publications. Users must ensure that the product(s) is appropriate for their application, prior to use. The information published in this publication is based on research and development work carried out in our laboratory and is to the best of our knowledge true and accurate. Micromaster Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are intended for laboratory, diagnostic, research or further manufacturing use only and not for human or animal or therapeutic use, unless otherwise specified. Statements included herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

