






## CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: Sterile Transport Swab	
Product Code: PW007-100PC	Expiry Date: 28/12/2021
Lot Number: 09071216	Release Date: 28/12/2016
Analytical Report Number: QAD/AR-16/12-1673	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Transparent, clear & unbreakable, polypropylene tube with specially designed cap and cotton bud with polypropylene stick. Individually packed, Size: 150mm x 12mm diameter. Pack of 100 nos.	Transparent, clear & unbreakable, polypropylene tube with specially designed cap and cotton bud with polypropylene stick. Individually packed, Size: 150mm x 12mm diameter. Pack of 100 nos.	Complies
Sterilization	The above batch of swabs was sterilized by Microtrol Sterilization Services Pvt. Ltd. by treatment with Ethylene oxide to the following parameters:  Concentration of Ethylene oxide : NLT 750 gm/M3 Exposure time : 6 hours Humidity : Min. 50% RH Temperature : 50 Deg C	Concentration of Ethylene oxide: 750 gm/M3	Complies
Sterility Testing	Test result should comply with sterility test (71) of USP 37 method.	Sterility passes	Complies

### Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
<b>TESTED BY</b> (Divya Sasidharan) MICROBIOLOGIST	<b>CHECKED BY</b> (Vijayalakshmi Poojari) HEAD MICROBIOLOGY	<b>APPROVED BY</b> (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report

Rev. 00

10/12/2014

QAD/FR/01

**Micromaster Laboratories Pvt. Ltd.**  
Manufacture & Exporter

Culturing your needs in Microbiology...

38/39, Kalpataru Industrial Estate, Opp Lawkim, Off. Ghodbunder Rd.,  
Near Runwal Estate, Thane (W) 400 607, Maharashtra, INDIA.  
Ph: +91-22-2589 5505, 4760, 4681, 9320126789  
Email: micromaster@micromasterlab.com

