






CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: Autoclavable Petri plates-100mm	
Product Code: PW008-10PC	Expiry Date: 05/12/2021
Lot Number: 08921216	Release Date: 05/12/2016
Analytical Report Number: QAD/AR-16/12-1593	Page No: 1 of 1

TEST CONDUCTED	SPECIFICATIONS	OBSERVATIONS	REMARKS
Appearance	Polypropylene, clear, translucent & Unbreakable.	Polypropylene, clear, translucent & Unbreakable.	Complies
Size	100mm diameter x 15mm	100mm diameter x 15mm	Complies
Use	Used for Microbiological testing.		

Final Results: **COMPLIES**

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
TESTED BY (Divya Sasidharan) MICROBIOLOGIST	CHECKED BY (Vijayalakshmi Poojari) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report

Rev. 00

10/12/2014

QAD/FR/01

Micromaster Laboratories Pvt. Ltd.
Manufacture & Exporter

Culturing your needs in Microbiology...

38/39, Kalpataru Industrial Estate, Opp Lawkim, Off. Ghodbunder Rd.,
Near Runwal Estate, Thane (W) 400 607, Maharashtra, INDIA.
Ph: +91-22-2589 5505, 4760, 4681, 9320126789
Email: micromaster@micromasterlab.com

