






CERTIFICATE OF ANALYSIS, QUALITY AND CONFORMITY

Product Name: Transport Swabs with Viral Transport Medium	
Product Code: TM006	Expiry Date: 31.07.2020
Lot Number: 00460119	Release Date: 18/02/2019
Analytical Report Number: QAD/AR-19/01-94	Page No: 1 of 1

PREPARED MEDIA			
TEST CONDUCTED	SPECIFICATION	OBSERVATIONS	REMARKS
Appearance	Sterile transport swab with Viral transport medium	Sterile transport swab with Viral transport medium	
Color of Final Medium	Light Pink colored Medium	Light Pink colored Medium	Complies
Clarity of Final Medium	Clear Solution	Clear Solution	Complies
pH of Final Medium at 25° C	7.4-7.8	7.49	Complies
Quantity of Medium	3ml	3ml	Complies
PREINCUBATION CHECK OF PLATE			
Pre Incubation Check	No growth should be observed after incubation at 30-35 °C for 24 - 48 hrs	No growth observed after incubation at 30-35 °C for 24 - 48 hrs	Complies

Final Results: COMPLIES

This is to certify that this lot passes and it conforms to the above mentioned tests and specifications.

		
TESTED BY (Pooja Sutar) MICROBIOLOGIST	CHECKED BY (Madhavi Sawant) HEAD MICROBIOLOGY	APPROVED BY (Rupali Sabale) HEAD QUALITY ASSURANCE

End Of The Report

TM006COA,QAD/FR/031,Rev.00/01.01.2018

Micromaster laboratories Pvt.Ltd
Manufacturer and Exporter
Unit No.38/39, Kalpataru Industrial Estate .Near Runwal Estate,Behind R Mall,
Ghodbunder Road,Thane (W)-400607,Maharashtra,India
[Tel:- +91-22-2589 5505](tel:+91-22-25895505)sales@micromasterlab.com